







# EnergyLight

HF3309, HF3308

**PHILIPS**

**ENGLISH 6**

**DEUTSCH 10**

**FRANÇAIS 15**

**NEDERLANDS 19**

**ITALIANO 23**

**NORSK 27**

**SVENSKA 31**

**SUOMI 35**

**DANSK 39**

## Introduction

The appliance is suitable for both light baths and for normal lighting. Light baths will help counteract the tiredness, low spirits and the lack of energy that many people experience in autumn and winter, when the days are darker and shorter. These symptoms, which are thought to be caused by a disruption of our biological clock, are often referred to as the winter blues. A shortage of light can also have a negative effect on our sleep-wake rhythm.

EnergyLight and the winter blues:

The winter blues are characterised by a combination of symptoms: tiredness, lack of energy and a need for more sleep, unexpected mood changes, having difficulty getting out of bed in the morning and having trouble concentrating. Sometimes the winter blues can also cause you to gain weight, as your body craves food that is rich in sugar and carbohydrates, such as pasta and chocolate. The winter blues occur during the darker period of the year, i.e. from September until March. More than 60% of the population is troubled by one or more symptoms of the winter blues. Many health experts agree that the winter blues can be counteracted by exposure to bright light. 'Bright light' means daylight with a certain intensity. In winter, the light intensity is often much lower than in summer: on a gloomy winter day the light intensity will vary between 2,500 and 10,000 lux, whereas this may be as much as 100,000 lux on a bright summer day. Moreover, in winter most people spend a lot of time indoors, either at home or at work, where artificial light often has an intensity of as little as 500 lux and sometimes even less. Together with scientists and medical specialists in the field of treating the winter blues, Philips has developed this special lighting appliance. This appliance allows you to safely get the amount of light you need to make the winter blues disappear. Clinical tests have shown that bright light works. In approx. 90% of all cases, the winter blues disappear when the EnergyLight is used regularly. Taking light baths is generally accepted as one of the safest and most effective ways of fighting the winter blues.

EnergyLight and the sleep-wake rhythm

A shortage of light plays a major role in the disruption of our sleep-wake rhythm. Hormones are no longer produced at exactly the right time of day, it becomes harder to sleep soundly throughout the night and the need for short naps during the day increases. Research has shown that an extra dose of light reinforces our sleep-wake rhythm and thus improves the quality of our sleep as well as the way we function during the day. The EnergyLight provides this extra dose of light.

For more information about the health benefits, see [www.philips.com](http://www.philips.com), search word 'EnergyLight'.

For more information about the use of the EnergyLight to counteract the negative effects of shift work, you can contact our specialised distributor:

MediluX BV, tel.: 0411-642834

## General description (fig. 1)

- A Energy lamps
- B Lamp housing with angle adjustment
- C User information
- D On/off button (HF3308 only)
- E Control panel (HF3309 only)

- 1 = on/off button
- 2 = light intensity ring
- 3 = timer
- 4 = energy level indicator
- 5 = distance selection ring
- 6 = pause button

### Important

Read these instructions for use carefully before using the appliance and save them for future reference.

- ▶ Check if the voltage indicated on the appliance corresponds to the local mains voltage before you connect the appliance.
- ▶ Keep this appliance away from water! Do not use this appliance in damp environments (e.g. in a bathroom or near a swimming pool).
- ▶ If the mains cord is damaged, it must be replaced by Philips, a service centre authorised by Philips or similarly qualified persons in order to avoid a hazard.
- ▶ Always unplug the appliance after use.
- ▶ Keep the appliance out of the reach of children.
- ▶ The duration of the light bath depends on the distance at which the appliance is used. See chapter 'Using the appliance' for the appropriate distance and duration. Repeat the light bath on at least five successive days (take them at a time that suits you and make this a daily habit during the course). You can start a new light bath course (5 or more successive days) whenever you feel you need it.
- ▶ Use the light of this appliance as an addition to the normal light (sunlight, electric light) in your home. Do not use the appliance in a room that is not illuminated by any other light source.
- ▶ It is not advisable to look into the light too long.
- ▶ After the first few light baths you may experience some eye strain and/or a light headache. These side effects tend to disappear quickly.
- ▶ Consult your doctor before you start using the EnergyLight if:
  - you have suffered from or are suffering from severe depression;
  - you are hypersensitive to light (e.g. diabetics and epileptics);
  - you have eye complaints;
  - you are on specific medication (i.e. certain antidepressants, psychotropic drugs or malaria tablets);
  - you have high blood pressure.
- ▶ If using the appliance causes you to become too active, we advise you to reduce the light dose. If there is no improvement, we advise you to consult your doctor.
- ▶ If you feel exhausted, distressed or restless and/or suffer from sleeping problems (e.g. insomnia) after having used the appliance, it is always advisable to consult your doctor or psychotherapist.
- ▶ This appliance meets the requirements of MDD 93/42/EEC for medical appliances.

### Preparing for use

- 1 Place the appliance on a tabletop.
- 2 Put the plug in a wall socket (220-240V) (fig. 2).
- 3 Switch the appliance on by pressing the on/off button.

**HF3308:**

The advisable duration of the light bath depends on the distance at which the appliance is used:

- approx. 2 hours at a distance of 60 cm
- approx. 1 hour at a distance of 40 cm
- approx. 30 minutes at a distance of 20 cm

**HF3309:**

- 1** Set the distance by turning the distance selection ring to the appropriate position.
  - 2** Set the light intensity by turning the light intensity ring to the appropriate position.
- When you set the distance selection ring to 20, 40 or 60 cm, the display shows the duration of the light bath. You can always adapt this time by setting a different light intensity. Lowering the light intensity leads to a longer duration and vice versa. For example, when you select the maximum light intensity and a distance of 20 cm, the duration of the light bath will be 30 minutes (fig. 3).
  - It is possible that the display does not show a light bath duration. This indicates that the set distance and light intensity will not yield the positive effects that can be expected if the prescribed combination of distance and time is observed. This may for instance occur if you select a too long distance (60 cm) in combination with a low light intensity.
  - When you set the distance ring to the manual setting, you can use the appliance when you are unable to sit at a set distance from it. At this setting you can adjust the light intensity but the display will not show a light bath duration. Because no distance has been selected, the effect cannot be measured or guaranteed.

When you switch the appliance on, it will start operating at the highest light intensity. After a few seconds, the light intensity will adjust itself to the level you have selected.

### Using the appliance

It is advisable to start using the EnergyLight as soon as you feel the winter blues coming on or when you feel like you need extra light.

- 1** Position yourself near the appliance and turn the screen in such a way that you can comfortably look in the direction of the light (fig. 4).

Make sure there is a distance between your face and the screen of 20 to 60 cm.

- You do not need to sit right in front of the appliance (fig. 5).
- **HF3309 only:** The energy level indicator shows the total amount of energy that has already been emitted during the light bath. The indicator is divided into ten parts that each represent 10% of the emitted energy.
- You do not have to take the whole light bath without interruption. If it suits your schedule better, you can interrupt your light bath and continue later.

- 2** **HF3309 only:** Press the pause button if you want to interrupt the light bath. The pause symbol in the display starts blinking. Press the pause button again if you want to continue with the light bath.

- 3** Look into the light every now and then.

Do not look into the light continuously! You can simply engage in other activities such as reading, writing, embroidering or handicraft while taking a light bath and look into the light every now and then.

- 4** Preferably take your light baths between 6 o'clock in the morning and 8 o'clock in the evening (fig. 6).

- ▶ *Note: if you want to readjust your biological clock, we advise you to take an extra light bath in the morning, immediately after waking up.*

**5 Repeat the light baths until you feel better.**

Generally, improvements will occur within 1 week.

You may repeat the light bath as often as you like during the dark period of the year. Take a light bath on at least 5 successive days.

**6 HF3308: Switch the appliance off after use and remove the mains plug from the wall socket.**

- ▶ **HF3309:** Approx. 2 minutes before the end of the light bath, the time on the display blinks during 1 minute. After this minute, the energy level indication also starts blinking until the appliance automatically switches off.

**7 Remove the mains plug from the wall socket.****Cleaning the appliance**

Always unplug the appliance and let it cool down before cleaning it.

**1 Clean the outside of the appliance with a damp cloth (fig. 7).****2 Make sure that no water runs into the appliance (fig. 8).**

Do not use abrasives, scourers or liquids such as alcohol, petrol or acetone.

**Replacement****Tube lamps**

- ▶ When the tube lamps stop working or after a period of 10,000 burning hours (after many years of normal use), they have to be replaced. Tube lamps must always be replaced by Philips or a service centre authorised by Philips, as special tools and/or parts are required.

The tube lamps must always be replaced with Philips Original EnergyLight PL-L 36W tube lamps.

**Environment**

Tube lamps contain substances that may pollute the environment. The lamps of the appliance contain substances that can be harmful to the environment. When you discard the lamps, do not throw them away with the normal household waste but hand them in at an official collection point.

- ▶ **Do not throw the appliance away with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this you will help to preserve the environment (fig. 9).**

**Guarantee & service**

If you need information or if you have a problem, please visit the Philips website at [www.philips.com](http://www.philips.com) or contact the Philips Customer Care Centre in your country (you will find its phone number in the worldwide guarantee leaflet). If there is no Customer Care Centre in your country, turn to your local Philips dealer or contact the Service Department of Philips Domestic Appliances and Personal Care BV.

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**PHILIPS DAP B.V.**