# beurer

BEURER GmbH + Co. KG • Söflinger Str. 218 • D-89077 Ulm (Germany)  
Tel. 0731 / 39 89-144 • Fax: 0731 / 39 89-145 • Mail: kd@beurer.de

## Diagnose-Waage BG 65


Gebrauchsanleitung  
Instructions for Use  
Mode d'emploi  
Instrucciones para el uso  
Istruzioni per l'uso  
Gebruikshandleiding  
Instruções de utilização  
Kullanma Talimatı  
Οδηγίες χρήσης

Diagnose-Waage  
Diagnostic scales  
Pèse-personne impédancemètre  
Báscula para diagnóstico  
Bilancia diagnostica  
Personenweegschaal  
Balança de diagnóstico  
Diyagnoz terazisi  
Διαγνωστική ζυγαριά

CE

# BG 65

## 🇩🇪 Deutsch

1. Display
2. „Auf“-Taste
3. „Ab“-Taste
4. Set-Taste
5. User-Taste
6. Elektroden

## 🇬🇧 English

1. Display
2. Up key
3. Down key
4. SET key
5. User key
6. Electrodes

## 🇫🇷 Français

1. Ecran
2. Touche Haut
3. Touche Bas
4. Touche SET
5. Touche User
6. Electrodes

## 🇪🇸 Español

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basierend auf: "Principles + Labs for Physical Fitness and Wellness. 1<sup>st</sup> edition by Copyright 1999"

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Dear customer,

thank you for choosing a product of our range. Our name stands for high and thoroughly checked quality-products in the field of warmth, gentle therapy, diagnosis/blood pressure, massage and air.

Please read the instruction for use carefully and pay attention to the instructions.

With regards

Your Beurer-Team

## Worth to know

### The measurement principle of the body water and fat scales

The body water and fat scales measures body water and fat by the principle of B.I.A., that is, bio-electric impedance analysis. An imperceptible and completely safe and harmless current measures the proportions of body tissues within seconds. Through this measurement of electrical resistance (impedance), the percentage of fat and also the percentage of water in the body can be determined. Muscle tissue and water have good electrical conductivity and thus low resistance. Fat tissue, in contrast, has low conductivity, as the fat cells hardly conduct the current because of their very high resistance.

### Transparent electrodes integrated in the platform

New and innovative technologies allow transparent glass coatings to be manufactured that are conductive at the same time. These act as electrodes for measuring body fat/body water.

The advantage of this technology is apparent in the design and cleaning is much easier compared to conventional electrodes.

### Weighing correctly

If possible, weigh yourself at the same time of day (ideally in the morning), after going to the toilet, fasting and without clothes, in order to obtain comparable results.

Only the long-term trend counts. Short-term fluctuations in weight within a few days are usually due only to loss of fluid. However, body water plays an important part in well-being.

The interpretation of the results depends on the changes in the

1. total weight
2. percentage of fat
3. percentage of water

and on the length of time over which these changes occur. Rapid changes over days should be distinguished from medium-term changes (over weeks) and long-term changes (months).

As a basic rule, SHORT-TERM weight changes represent changes in water content almost exclusively, while medium- and long-term changes can ALSO refer to the fat content.

- If the weight drops in the short-term but the proportion of body fat increases or remains the same, you have only lost water, e.g. after a training session, going to the sauna or after a diet limited only to rapid weight loss.
- If the weight increases in the medium term, while the proportion of body fat falls or remains the same, you might have built up your muscle mass or there can be increased retention of body water (e.g. as a result of taking medications or cardiac function becoming weaker in elderly people)
- If weight and body fat proportion fall at the same time, your diet is working – you are losing fat mass.

Medium-term changes in the proportion of body water can also be due to changes in muscle mass.

# ! Important Notes – keep for later use!

Warning: The scales must not be used by persons with medical implants (particularly heart pacemakers)! Otherwise their function may be impaired.

## 1. Caution – For your safety

- Do not use on pregnant woman. Inaccurate measurements may arise due to the amniotic fluid.
- Caution! Don't step on the scales with wet feet and don't step on the damp surface – Slippery surface!
- Protect your personal scales from impact with hard objects, moisture, dust, chemicals, toiletries, liquid cosmetics, great temperature fluctuation and closeness to sources of heat (open fires, radiators).
- Repairs may only be performed by Beurer Customer Service or by accredited retailers. However, always test the batteries and replace them if necessary, before returning for repairs. If repair work is necessary, however, first remove the batteries and then pack the personal scales in the original box prior to transportation.
- All our personal scales comply with EU guideline 89/336 + relevant supplements. Should you have any questions regarding operation and usage of our appliances, please contact your local retailers or Beurer Customer Service.
- Clean the housing of your personal scales with a soft and slightly damp cloth. Do not use any aggressive detergents and do not submerge the scales.
- All personal scales are intended for private use and have not been developed for medical or commercial purposes.
- An age range of 10 ... 100 years and heights of 100 ... 220 cm (3-03" – 7-03") can be preset for the body fat analysis. The scales can accept a maximum load of 150 kg (330 lb, 24 St). The weight measurement is displayed in 100 g (0,2 lb, 1/4St) divisions. The measuring results for body fat are displayed in 0.1% divisions. The results of the body fat and water analysis are displayed in 0.1% steps.
- On delivery, the scales are set to the "cm" and "kg" units. The switch at the back of the scales allows you to change to "inch" (") and "pound/Stones" (lb, St).

## 2. Operation

### 2.1 Batteries


Pull off the battery insulating strip on the cover of the battery compartment and remove the battery's protective film and insert the battery according to polarity. If the scales show no function remove the batteries completely and put them in again.

Your personal scales are fitted with a "replace battery" display. If the personal scales are operated using an extremely weak battery, the message "Lo" will appear on the display panel and the personal scales switch off automatically. In this case, the batteries must be replaced (1 piece 3V lithium battery CR 2032).

Spent batteries and rechargeable batteries do not constitute normal household waste! They are considered to be toxic waste and, as such, should be disposed of in special containers, toxic waste collection points or brought to electrical goods dealers. You are legally obliged to dispose of the batteries.

Note: Batteries containing pollutant substances are marked as follows: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

Position the personal scales on a firm, flat surface; a firm flooring is a basic requirement for correct weight measurement.

## 2.2 Measuring weight only

Tap the platform of your scales briefly but firmly with your foot.

Briefly touch the scale with your foot. The whole display (fig. 1) appears as an automatic test until „0.0“ (fig. 2) is displayed (tap-on technology). The scale is now ready to measure your weight. Stand on the scale without moving about and distribute your weight evenly on both legs. The scale immediately starts to measure your weight.

„000“ flashes until your weight has been measured (fig. 3).

Your measured weight is displayed shortly thereafter (fig. 4). The scale switches itself off a few seconds after you step down from it.

fig. 1

fig. 2


fig. 3

fig. 4


fig. 5

fig. 6

fig. 7

fig. 8


fig. 9


fig. 10

## 2.3 Measuring weight and body fat and body water percentage

Your personal parameters must first be stored.

– With BG 65 you have to step briefly onto the platform (2.2) and then press "SET". The memory places 1 to 10 now appear in the display (fig.5). Select the desired memory place by operating the "up" ▲ or "down" ▼ buttons and confirm the desired memory place with the "SET" button.

– Height („Ht“) in „cm“ (inch) (fig. 6) now blinks in the display.

Select the correct height by either repeatedly or continuously pressing „up“ ▲ or „down“ ▼ and confirm by pressing „SET“.

– „Age“ in years flashes in the display (fig. 7).

Select the correct age by either repeatedly or continuously pressing „up“ ▲ or „down“ ▼ and confirm by pressing „SET“.

– The sex setting „Woman“ or „Man“ appears on the display (fig. 8). Select the sex by operating the "up" ▲ or "down" ▼ buttons and confirm the desired memory place with "SET".

– The setting for your individual degree of activity and fitness now appears on the display (fig. 9–11). Select the desired level of activity by operating the "up" ▲ or "down" ▼ buttons and confirm the desired memory place with the "SET" button.

**Fitness 1:** No or little physical activity

(less than 20 minutes of physical exercise twice a week, e.g. going for a walk, light gardening, gymnastic exercises)

**Fitness 2:** Moderate physical activity

(physical activity up to 2 – 5 times a week, 20 minutes each time, e.g. jogging, cycling, tennis, ...)

**Fitness 3:** Intensive physical activity

(Daily intensive training or physical work, intensive running, construction worker, ...)

– The values, which have now been set, are displayed once again, one after the other. The scale then automatically turns itself off.

The scale has 10 memory spaces for basic data, where personal details, of individual family members for example, can be stored and called up at any time.

– You have to switch on the scales first by stepping on the platform. Select the memory place where your basic personal data are stored by pressing the "User" button several times. The data are displayed in succession until "0.0" appears (fig. 2).

– Step on the scales barefoot and ensure that you are standing on the two electrodes. Your body weight is first determined and displayed as described under 2.2.

– Remain standing on the scales and the body fat and body water analysis is now performed. This can take a few seconds. During this time "bF" and "H2O" flash alternately on the display (fig. 12 + 13).

**Important:** Your feet, legs calves and thighs must not touch each other, as this may result in your body fat being measured incorrectly.

- Finally your body fat proportion ("bf"). and your body water proportion ("H2O") are displayed in %. The scales switches off automatically after approx. 5 seconds.


fig. 11

## 2.4 Tips for using the scale

Important when measuring body fat/body water:

- The measurement must only be made barefoot and it is helpful if the soles of the feet are slightly damp. Completely dry soles can lead to unsatisfactory results, as these have too little conductivity.
- Stand still during the measurement
- Wait a few hours (6-8) after unusually strenuous activity
- Wait for about 15 minutes after getting out of bed so that the water in your body can become distributed.

fig. 12

fig. 13

Measurement of body fat and body water is not reliable in:

- children under 10 years
- individuals with fever, undergoing dialysis, with symptoms of edema or osteoporosis
- individuals taking cardiovascular medication. Persons taking vasodilating or vasoconstricting medications
- Persons with substantial anatomical deviations in the legs relative to total height (leg length considerably short or long).

## Body fat guide

The following body fat levels provide you with a guideline (for further information, please consult your doctor.) A lower level is often found in athletes. Depending on the type of sport, the intensity of training and physical constitution, levels can be obtained that are even lower than the stated guidelines.

### Women

Age	very good	good	average	poor
<19	<17%	17-22%	22,1-27%	>27,1%
20-29	<18%	18-23%	23,1-28%	>28,1%
30-39	<19%	19-24%	24,1-29%	>29,1%
40-49	<20%	20-25%	25,1-30%	>30,1%
>50	<21%	21-26%	26,1-31%	>31,1%

### Men

Age	very good	good	average	poor
<19	<12%	12-17%	17,1-22%	>22,1%
20-29	<13%	13-18%	18,1-23%	>23,1%
30-39	<14%	14-19%	19,1-24%	>24,1%
40-49	<15%	15-20%	20,1-25%	>25,1%
>50	<16%	16-21%	21,1-26%	>26,1%

based on: "Principles + Labs for Physical Fitness and Wellness. 1<sup>st</sup> edition by Copyright 1999"

According to the World Health Organisation (WHO 2001), the percentage of body water should lie within the following ranges:

Women: 50 – 55 %

Men: 60 – 65 %

Children: 65 – 75 %

In persons with a high proportion of fat, the percentage of body water obtained is very often below the stated guidelines.

### 3. Incorrect measurement

If the scale detects an error during weighing, „Err“ appears in the display.

Possible causes for errors:	Error handling:
- The scale was not activated before being stepped on to	- Activate the scale correctly and reweigh
- The maximum 150 kg weight-bearing capacity was exceeded	- Only weigh maximum permitted weight
- The electrical resistance between the electrodes and the soles of your feet is too high (e.g. in cases of extremely hard skin)	- Please repeat the measurement with bare feet. If necessary, you can wet your feet a little. If necessary, calluses should be removed on the soles of the feet.
- You did not stand still enough on the scale	- Please stand still
- Your body fat is not within the measurable range (less than 5% or greater than 50%)	- Please weigh again, this time with bare feet. If necessary, you can wet your feet a little.
- The percentage of water lies outside the measurable range (less than 43% or greater than 75%)	- Please weigh again, this time with bare feet. If necessary, you can wet your feet a little.

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Baseado em: "Principles + Labs for Physical Fitness and Wellness. 1ª edition by Copyright 1999"